

Ultrasonic Cleaners Cleaning Solutions & Tablets



CLEANLETS - CLEAN TABLETS - SULTAN

Features:

Every Cleanlets ultrasonic tablet provides quick, effective cleaning in a single dose. Patented technology gives you highly concentrated cleaners that are safe and effortless to use. Effervescent formula makes mixing easy and automatic. Its unique blister packaging makes Cleanlets a snap to dispense.



Code	Name	Price	+
402344	SULTAN - Cleanlets Ultrasonic Instrument Cleaning Tartar & Stain Tablets - 32/b...	US\$ 68.39	+Info
402343	SULTAN - CleanLets Ultrasonic Cleaning Tablets- 32/bx #SUL 21500	US\$ 68.39	+Info

VISION - 121460 Hyper Clean For Metals Gallon

US\$ 25.81

VISION - Hyper Clean For Metals Gallon - # Hcm [VISION - Hyper Clean Metal Loupes & Magnifiers](#)
121460 - Mfg: HCM

Descriptions

[HyperClean for Metal.](#)

HyperClean Metal is NEW Technology. Ultrasonic solution for the dental laboratory that completely eliminates contamination on Metal / Dental Alloy surface.

NASA technology

- * Completely eliminates contaminants
- * Helps eliminate discoloration
- * Removes radioactivity
- * Bio-degradable
- * Non-toxic
- * Non-carcinogenic
- * Avoiding Risk in Metal restorations

It has been said that the ability to avoid risk, is a great insurance policy in any marketplace. This is especially important in a doctor/patient relationship when restorative dentistry is part of the equation.

Reducing risk can pay large premiums when utilizing new technology that has the possibility to dominate a market, but is unproven over time.

The Zirconia/ceramic restoration is such a technology. The esthetic and health benefits are becoming well accepted, as well as the stable pricing



from the laboratory.

Quantity: 1 gallon

**Cory Labs - 400672 ProOne Sonic Cleaner
tabl 80pk**

US\$ 44.93

Cory Labs - ProOne Sonic Cleaner Tablets 80/bx.
#POS80Cory Labs - ProOne Sonic Cleaner Tablets
Infection Control
400672 - Mfg # POS80

Features:

Pro One Sonic Tablets Ultrasonic Cleaner by Cory Labs is the only bio-enzymatic tablet on the market. It will degrade all complex organics and soils that traditional protease enzymatic formulas cannot.

- * Ultrasonic Tabs by Cory Labs are easy-to-use enzymatic tablets.
- * Dissolve in a gallon of water in under a minute to create a powerful and effective ultrasonic cleaning solution.
- * Made from a biodegradable, non-chlorine, non-corrosive formulation that contains fast acting protease enzymes.
- * "The solution will greatly assist in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents.
- * bold; "Virtually eliminates handling and storage of heavy liquid concentrates.
- * Contains a pleasant fragrance and is a neutral pH. May also be used as an instrument pre-soak or as an evacuation system cleaner.

Quantity: 80 Tablets per box



Cory Labs - 400673 Purity Water Line Cleaner 60pk

US\$ 22.43

Cory Labs - Purity Water Line Cleaner 60/bx #PRT60Cory Labs - Purity Water
Line Cleaner
cleaning Solution & Tablets
400673 - Mfg # PRT60

Features:

Purity Waterline Cleaner is a safe and effective procedure to maintain clean dental unit waterlines.

- * Contains natural preservatives for continuous cleaning of self-contained water systems



- * **Environmental friendly, no harsh chemicals.**
- * **Effective, easy to use, no purging required**
- * **Orderless, tasteless, non-irritating and has no effect on bonding**
- * **Compatible with all dental unit water systems.**

Quantity: 60 Tablets per box

Cory Labs - 400674 Stain Out Tartar & St Remover 40pk

US\$ 27.38

Cory Labs - Stain Out Tartar & Stain Remover Tablets 40/bx #TSR40
Cory Labs - Stain Out Tartar & Stain Remover Cleaning Solution & Tablets
400674 - Mfg # TSR40



Features:

Stain Out - Tartar and Stain Remover Tablets 40/bx - Cory Labs are Fast and Easy Tablet uses high level formula is excellent for cleaning and removal of tarter, stain, tobacco, calculus and all other food stains. Can be used on dentures, bridges and all orthodontic appliances.

- * **Fast and Easy Tablet uses high level formula**
- * **excellent for cleaning and removal of tarter, stain, tobacco, calculus and all other food stains.**
- * **Can be used on dentures, bridges and all orthodontic appliances.**
- * **Can also be used for removal of permanent cements**

Quantity: 40 Tablets per box

SYMPRO CLEANING FLUIDS - RENFERT

SYMPRO CLEANING FLUIDS - RENFERT

SYMPROfluid Universal

Is used for cleaning in combination with the device. This fluid is used for removing "standard" accretion such as plaque, calculus, discoloration and adhesive denture creams.

SYMPROfluid Nicoclean

Also removes heavy nicotine stains from dentures. It is used in the form of follow-up treatment without the use of the device.



Code	Name	Price	+
23-65000600	RENFERT - Sympro Fluid Universal 2L x 2- # 6500-0600 # 65000600 This product ca...	US\$ 58.75	+Info
23-65000610	RENFERT - Sympro Fluid Nicoclean 2 x 2L- # 6500-0610 # 65000610 This product ca...	US\$ 58.75	+Info

Ultrasonic Cleaners

Ultrasonic Accessories



ACCESSORIES FOR EASY CLEAN ULTRASONIC UNIT S - RENFERT

ACCESSORIES FOR EASY CLEAN ULTRASONIC UNITS - RENFERT

These accessories will make your portable dental units more efficient and improve the quality and function of your machine.



Code	Name	Price	+
120759	RENFERT - Front Foil For Easyclean - # 900030775	US\$ 69.53	+Info
23-1850-0005	RENFERT - Plastic Acid Bath Insert for EasyClean- # 1850-0005 # 18500005	US\$ 78.88	+Info
23-1850-0003	RENFERT - Stainless Steel Basket for EasyClean- # 1850-0003 # 18500003	US\$ 72.82	+Info
23-1850-0002	RENFERT - Stainless Steel Cover for EasyClean # 1850-0002 # 18500002	US\$ 66.73	+Info
23-1850-0004	RENFERT - Stainless Steel Immersion Basket for EasyClean- # 1850-0004 # 18500004	US\$ 66.73	+Info
23-1850-0001	RENFERT - Plastic Lid for EasyClean-# 1850-0001 # 1850-0001 # 18500001	US\$ 37.07	+Info
23-1850-0006	RENFERT - Cleaning Jar with Lid and rubber band 600ml for EasyClean-# 1850-0006 ...	US\$ 25.17	+Info
120758	RENFERT - Knob For Easyclean - # 900030774	US\$ 20.31	+Info
120757	RENFERT - Fuses For Easyclean - # 900030118	US\$ 15.38	+Info
23-1850-0007	RENFERT - Plastic container with Lid for EasyClean- # 1850-0007 # 18500007	US\$ 16.00	+Info
120701	RENFERT - Plastic Container For Easyclean - # 18500007	US\$ 17.56	+Info

RENFERT - 23-65000450 SYMPRO disinfectable beaker

US\$ 25.54

RENFERT -SYMPRO disinfectable beaker - #65000450RENFERT - SYMPRO disinfectable beaker

Solutions

23- 65000450 Mfg # 65000450



Features:

NEW! Operating and hygiene preparation guidelines for using the SYMPRO, [download here](#)

Technical data:

Power supply 100 - 240V / 50/60Hz

Power consumption 90VA
Dimensions (W x H x D) 150 x 240 x 280 mm
Weight w/o beaker ? 3kg
Noise level
Rotational speed 1/min

Scope of Delivery
SYMPRO cleaning beaker disinfectable

Ultrasonic Cleaners

Ultrasonic Units



EASYCLEAN ULTRASONIC CLEANER - RENFERT
EASY CLEAN ULTRASONIC CLEANER - RENFERT

Ultrasonic cleaner used for time saving and effective cleaning.
 Ultrasonic cleaning agent is the perfect partner for optimum cleaning.
 Ultrasonics is the most up-to-date and effective process available for fine cleaning.
 The Easyclean provides state-of-the-art ultrasonic cleaning technology for all types of ultrasonic cleaning.



Advantages

- * Varied range of applications.
- * Cleans items thoroughly and gently, even highly intricate or delicate items.
- * Saves time by eliminating manual cleaning.
- * Extremely easy to operate.

US DENTAL DEPOT is an Authorized RENFERT Dealer in US.

Code	Name	Price	+
23-1850-1000	RENFERT - EasyClean Ultrasonic Cleaner 110 Volts-# 1850-1000 # 18501000 Tank- in...	US\$ 855.00	+Info
23-1850-0000	RENFERT - EasyClean Ultrasonic Cleaner 220 Volts-# 1850-0000 # 18500000 Tank- in...	US\$ 879.36	+Info

ECCO - 112242 Ultrasonic Cleaner 0.6 Lt 110v

ECCO Ultrasonic Cleaner Stainless Steel Tank 600ml
Ultrasonic Units
112242

US\$ 63.00



The function of 600ml Ecco Ultrasonic cleaner
 The Ultrasonic cleaner machine uses water, alkaline or acidic water-based solvents as cleaning agents. Mainly applicable to clean plating parts before and plating after cleaning.

Features:

- * Stainless Steel Tank
- * moisture-proof PCB board
- * Quiet- Solid State Circuit
- * 3 min Auto Shut Off Timer
- * Temperature tunable from 20 to 80

Specifications:

- * Ultrasonic Power: 35Watts
- * Ultrasonic Frequency: 42000Hz
- * AC100-120V,60Hz
- * Tank size: 155 x 95 x 52mm
- * Overall size: 200 x 140 x 125mm
- * Volume: 600ml



ECCO - 101944 ECCO US-2.5 lit Digital Heated

US\$ 115.50

**ECCO - US-2.5 Digital Ultrasonic Cleaner
Ultrasonic Units
101944**

Features:

- * With Heater included
- * Large Industrial Grade Digital Ultrasonic Cleaner has powerful transducer, heater and large tank capacity for superior cleaning results, making it great for professional cleaning.
- * It is the first utilization of digital program in ultrasonic cleaner to control the unit, according to the quantity and condition of the items to be cleaner to choose a working cycle.

Descriptions:

- * High Power Transducer for Superior Cleaning Results
- * Isolated-Circuit Design, Keep Water Away for Enhanced Safety
- * First Ultrasonic Cleaner to Introduce Cooling Fan for Prolonged Operation

Applications

- * Great Cleaning Device for Scientific Labs, Medical and Dental Clinics, Electronics Workshops Hobbyists
- * Electronics Components, Mechanical Parts, Firearms,



Golf Club

* Personal Care Jewelry, Watches Straps and Bands, Glasses, Dentures, Toothbrushes, Table-ware

Specifications:

- * Ultrasonic Frequency: 42,000 Hz
- * Tank Material: Stainless Steel SUS304
- * Tank Capacity: 2500ml
- * Quiet, Solid States Circuitry
- * Timer: 5 Cycles, with Digital Time Display
- * 5 cycle timer : 1.5 min 3 min - 5 min - 6 min
- * Power Supply :AC 100 ~ 120V, 60Hz / AC 220 ~ 240V, 50 Hz
- * Power:160 W (110 V) / 170 W (220 V)
- * Unit Size: 260 x 270 x 140 mm (L x W x H) 11 x 9 x 6 inch
- * Inner Size: 250 x 160 x 85 mm (L x W x H) 9.84 x 5.90 x 3.15 inch
- * Box size: 340 x 255 x 225 mm (L x W x H) 13.39 x 10.04 x 8.86 inch
- * N.W. 2.5 kg
- * Large Tank: 2.5 Liters - 5.29pints - 160W - 110v

Shipping weight: 3 kg

SYMPRO DENTURE CLEANER - RENFERT

SYMPRO DENTURE CLEANER - RENFERT

Function & Performance

SYMPRO is a compact, high-performance cleaning unit for dentures, orthodontic appliances and splints.

Advantages

- * Maximum oral hygiene increases the well-being of the patient and promotes patient loyalty
- * The fully automatic process involves virtually no labor costs.
- * Time-saving work preparation for repairs, extensions and relines.



Code	Name	Price	+
23-65000000	RENFERT - Sympro Denture Cleaner with cleaning fluid- # 65000050 Rated voltage: ...	US\$ 877.55	+Info

ACCESSORIES

Code	Name	Price	+
23-65000550	RENFERT - Sympro Cleaning Needle 75gr Refill- # 6500-0550 # 65000550	US\$ 49.29	+Info
120781	RENFERT - Retainer Magnet (sympro) - # 900034039	US\$ 29.71	+Info
120780	RENFERT - Plastic Tweezers (sympro) - # 900034029	US\$ 18.28	+Info

ECCO - 100596 Ultrasonic Cleaner 1.4Lts 110V

US\$ 185.38

ECCO - Digital Ultrasonic Cleaner - 1.4lts/ 0.5gal - Stainless steel - 110volts- Heater included - Timer - 50wats consumption - Digital control and lecture - Stainless steel - [J#us71] 40kHz transducers

ECCO - Digital Ultrasonic Cleaner
Ultrasonic Units
100596



Features:

- * 1.4 lts / 0.5 gal
- * Heater included
- * Timer
- * 50 wats consumption
- * Digital control and lecture
- * 40 kHz transducers

Volts: 110v

Material: Stainless steel

ECCO - 101474 Ultrasonic Cleaner 3.5Lts 110V

US\$ 294.00

ECCO - Digital Ultrasonic Cleaner - 3.5 liters /0.80gal - Stainless steel - - 110volts- Heater included - Timer - 50wats consumption - Digital control and lecture - Stainless steel - [J# us32] 40kHz transducers- basket size - 10 1/2 in x 6 1/2in x 4in - instruments lay down in basket

ECCO - Digital Ultrasonic Cleaner
Ultrasonic Units
101474



Features:

- * 3.5 liters / 0.80 gal
- * Heater included
- * Timer
- * 50 wats consumption
- * Digital control and lecture
- * 40 kHz transducers
- * Basket size - 10 1/2 in x 6 1/2 in x 4 in
- * Instruments lay down in basket

Material: Stainless steel

Volts: 110v

ECCO - 100513 Ultrasonic Cleaner 600ml 110V

US\$ 115.50

ECCO - Digital Ultrasonic Cleaner -600ml - 0.7lts/ 0.25gal Stainless steel - 110volts- Heater included - 50wats consumption - Digital control - (J#US20 # PS-10A) [25kHz transducers]
ECCO - Digital Ultrasonic Cleaner Ultrasonic Units # 100513



Features:

- * 0.7 lts / 0.25 gal
- * Heater included
- * 50 wats consumption
- * Digital control
- * 25kHz transducers

Volts: 110v

Material: Stainless steel

ECCO - 114964 Ultrasonic Cleaner 7Lts 110V

US\$ 504.00

ECCO - Digital Ultrasonic Cleaner - 7 liters /1.60gal Stainless steel - - 110volts- Heater included - Timer - 50wats consumption - Digital control and lecture - Stainless steel - [J# US83) 40kHz transducers- - instruments lay down in basket -
ECCO - Digital Ultrasonic Cleaner Ultrasonic Units # 114964



Features:

- * 7 liters /1.60gal
- * Heater included
- * Timer
- * 50 wats consumption
- * Digital control and lecture
- * 40 kHz transducers
- * Instruments lay down in basket

Material: Stainless steel

Volts: 110v



Catalog Index / Indice

Cleaning Solutions & Tablets	1
Ultrasonic Accessories	4
Ultrasonic Units	5